

TESI Massage- and Traction System and accessories



TESI ComfoTrac duo Pt.-No. 018-0-0030



Basic functions

- Adjustable and regulated surface heating of the cover in the region of the back
- Comfortable 75 cm wide lying surface
- Vertebral massage with > adjustable massage area > adjustable pressure > adjustable speed
- Massage vibration with adjustable BMS frequencies
- All settings via a large LCD touchscreen
- Electrically adjustable stepped bed, subject to size of person
- Simple handling via the touchscreen display

Functions

- Cervical spine extension with precision-controlled traction force up to 160 N
- Lumbar spine extension with precision-controlled traction force up to 600 N

Indications

- Back pain
- Spinal disc herniation

Technical Data

Massage:	
Vibration frequency	15 – 35 Hz
Heating temperature	max 35°C
Treatment timer	1 – 30 min
Extension:	
Max traction force Lumbar spine	Continuous and intermittent max. 600 N
Cervical spine	max. 160 N
Power supply	100 – 240 V, 50 – 60 Hz
Class of equipment	IIb
Class of protection	I
Degree of protection	B
Dimensions	200 x 94 x 66 cm (W x D x H) height with arm 175 cm
Weight (without equipment)	140 kg
Warranty period	24 month

Subject to technical changes without prior notice.

Accessories

	Part-No.	Specification	Use	Application	Remarks
	018-0-0001	Thorax belt	Fixing the body while doing a traction on lumbar spine	On the thorax and loops at the head end	With TESI ComfoTrac and ComfoTrac Duo
	018-0-0002	Pelvic belt, Size M (90 – 110cm)	For lumbar spine extension	Around the waist and at the lumbar spine gallows	With TESI ComfoTrac and ComfoTrac Duo
	018-0-0003	Pelvic belt, Size S (75 – 90 cm)	For lumbar spine extension	Around the waist and at the lumbar spine gallows	With TESI ComfoTrac and ComfoTrac Duo
	018-0-0056	Pelvic belt, Size L (110 – 145 cm)	For lumbar spine extension	Around the waist and at the lumbar spine gallows	With TESI ComfoTrac and ComfoTrac Duo
	018-0-0040	Cervical carriage	For cervical spine extension	At the head and cervical gallows	With TESI ComfoTrac Duo
	018-0-0004	Cervical belt	For cervical spine extension (Alternative to cervical carriage)	At the head and cervical gallows	With TESI ComfoTrac Duo
	016-0-0058	Cervical cushion, small	For excessive curvature of the upper dorsal region	Is placed on the lying surface	With TESI ComfoTrac and ComfoTrac Duo

