

# NIDO

## BABY COMFORT PROGRAM

GINEVRI  
**70**  
1954-2024  
  
MADE IN ITALY



**GINEVRI**  
neo-tech for life

# Designed specifically for the comfort and the safety of the newborn.

The “Nido” baby crib is a modern solution designed with healthcare professionals in mind.

The accompanying trolley is strong and safe, made from durable epoxy-painted steel. This ensures both longevity and safety. The crib’s clever design allows the mother to stay close to her child, following the principles of **Family-Centered Care** and letting her stay comfortably in bed.

Safety is a priority with the rounded hood, providing protection for the newborn and the nurse. **Easy cleaning is guaranteed** with non-toxic materials, and the bassinet is versatile, easily adjusting to Trendelenburg and Fowler positions. The innovative handle on the crib’s hood allows for a full rotation, ensuring easy access to the newborn.

The crib’s transparent plastic structure not only makes cleaning simple but also allows for **phototherapy without moving the baby**. The forward-facing handle promotes easy movement and continuous observation during transport.

The standard “Nido” package includes a trolley with four castors (two lockable), a fire-retardant mattress, and an easily cleanable plastic bassinet.

More than just a crib, “Nido” transforms into a complete neonatal care system with various accessories. The exclusive Round Cover Hood and Screen accessories provide advanced comfort and **optimal protection during internal transport**.

Experience the next level of neonatal care with “Nido” – where safety, functionality, and Family-Centered Care come together for the well-being of both newborns and mothers.

## TECHNICAL SPECIFICATION

<b>Tilting (Trendelenburg and Fowler)</b>	+/- 12°	<b>hood</b>	3 kg
<b>Max. Load</b>	10 Kg	<b>Castors:</b>	4
<b>Dimensions</b>	80 x 50 x 92 cm	<b>Diameter:</b>	10 cm
<b>WEIGHT</b>		<b>Number of breaks:</b>	2
<b>trolley + cradle</b>	12 kg		

COVER SYSTEM

HANDLE

BABY FLAG  
(Optional)

BASKET (Optional)

WHEELS



## ACCESSORIES

Baby Bundle


Screen

Nido equipped with aquatherm,  
oxygen cylinder trolley, screen



Round Cover hood

Basket



# **G/NEVRI**

neo-tech for life

GINEVRI srl  
Via Cancelliera, 25/b  
00041 Albano Laziale (Rome) - Italy  
Phone +39 06 9345 100  
e-mail: mail@ginevri.com

Export Department  
Phone +39 06 93459 331  
e-mail: export1@ginevri.com

[ginevri.com](http://ginevri.com)

Authorized Distributor

The specifications in this catalogue are indicative. The company GINEVRI reserves the right to make changes, without further notice, to the products described within this catalogue in order to improve reliability, function or design.

Quality System  
ISO 9001:2015  
ISO 13485:2016

